23-315-9057

Nov 11, 2023 RECEIVED DATE Nov 09, 2023 SEND TO **35432** 



PAGE 1/2

| SSUE DATE |
| Nov 11, 2023

13611 B Street • Omaha, Nebraska 68144-3693 • (402) 334-7770 www.midwestlabs.com

# REPORT OF ANALYSIS

For: (35432) OPEN FARM Dog-OF Rawmix Fish

	Level Found		Reporting		Analyst-	Verified-
Analysis	As Received	Units	Limit	Method	Date	Date
Sample ID: 31023FWODRB	Lab Number: <b>14178709</b>	Date Sample	d: <b>2023-11-</b>	06		
Salmonella	negative	org/375g	1	AOAC OMA #2013.02; AFNOR QUA 18/08-03/15	smb8-2023/11/10	jzh4-2023/11/10
E. coli (generic)	n.d.	cfu/g	10	AOAC OMA 2018.13	Jhv0-2023/11/10	snl7-2023/11/10
Listeria monocytogenes	negative	org/125g	1	Real-Time PCR Assay; AOAC-RI 121402/AFNOR QUA 18/10-01/19	tma2-2023/11/11	jzh4-2023/11/11

For questions please contact:

Jamie Wood Account Manager

jwood@midwestlabs.com (402)590-2964

23-315-9057

Nov 11, 2023
RECEIVED DATE
Nov 09, 2023

SEND TO **35432** 



PAGE 2/2

| SSUE DATE |
| Nov 11, 2023

13611 B Street • Omaha, Nebraska 68144-3693 • (402) 334-7770 www.midwestlabs.com

## REPORT OF ANALYSIS

For: (35432) OPEN FARM Dog-OF Rawmix Fish

# **Detailed Method Description(s)**

## E. coli and Total Coliform using 3M Petrifilm

Sample analysis follows MWL MI 292 which is based on AOAC 2018.13. A representative sample is obtained and added to phosphate buffer. Aliquots of the sample are withdrawn and placed on Petrifilm plates. The plates are incubated for 18 to 24 hours. After incubation, the plates are counted to determine the number of generic E. coli and total coliforms present. The color of the colony and the presence of gas differentiate a generic coliform from E. coli. The levels are reported as colony forming units (cfu).

### L. mono real-time PCR

Sample analysis follows MWL MI 343 which utilizes the BAX® System Real-Time PCR Assay for Listeria monocytogenes. A representative sample is added to selective media, then incubated. After incubation, the sample DNA is analyzed using the BAX® System. Results are reported as negative or presumptive positive for L. monocytogenes.

### Sal - real time PCR

Sample analysis follows MWL MI 334 which is based on AOAC OMA 2013.02. A representative sample is obtained and combined with enrichment media, then incubated. After incubation, DNA is extracted from the sample and analyzed using BAX Real-Time PCR. Results are reported as negative or presumptive positive for Salmonella.